Preparation, Specimen Collection, Test Procedure, and Interpretation of Results

The collector should be trained on the procedures for operating the test. Do not break the seal of the protective pouch until ready to begin testing. Visually inspect the foil package to ensure it has not been compromised before beginning the test. If the package is not intact, the integrity of the test cup or dip may be compromised.

A fresh urine specimen should be collected directly into the cup and does not require any special handling or pre-treatment.

Test cups employ a temperature strip to validate that the urine is a freshly-voided specimen and should be checked immediately after collection.

NOTE: Urine specimens can be transferred from a urine collection container into the test cup, if necessary.

Procedure

1. Have a Test Results Record form ready to record results.
2. Require the donor to present identification.
3. Ask donor to remove any unnecessary outer clothing and empty all pockets.
4. Keep all backpacks and purses out of the restroom. Insure that all backpacks, wallets and purses are secured, preferably in a locked cabinet. If using a lockable cabinet or container, make sure that the donor is given the key.
5. Secure the collection site (restroom). Turn off water sources (sink faucets), ensure blue toilet water, inspect site for and remove cleaners and other substances, secure toilet tank top or blue tank water.

6. Tear open the foil pouch and remove the test.

7. Issue a cup to the individual to be tested. Instruct the donor not to flush the toilet or turn the water on and to open the door when they are done.

8. Have donor void urine directly into the test cup or specimen collection cup. Ensure that the specimen has a volume of at least 30 ml.

9. Check temperature (normal range of 90° to 100°F or 32° to 38°C) and observe the specimen for foreign material and/or discoloration.

10a. DIP TEST: Remove cap. Dip test strips into specimen for 15 seconds. Replace cap.

10b. TEST CUP: Pull the privacy label off to view results.


11. Read drug test strip results at 5 minutes, adulterants at 1 minute, and alcohol strips at 2 minutes, if applicable. Have donor sign and date the label on the test device. To interpret specimen validity and alcohol tests, ensure that the colour reference chart is present to properly analyze the test results for each of the test parameters.

12. Fill out Test Results Record with complete donor and test information.


13. If laboratory-based confirmation test is required, for non-negative result, secure specimen in appropriate shipping container and send to a confirmation laboratory with a signed Custody and Control Form (CCF).

Drug Test Results


Non-negative: One coloured line only visible in each Control Zone. No coloured line whatsoever appears in the Test Zone, indicating a Non-negative result for the corresponding drug of that specific Test Zone. It is recommended that the urine specimen is sent to a certified laboratory for confirmation testing.

Negative: Two coloured lines visible in each Control Zone and in the Test Zone, indicating that the concentration of the drug is below the detection limit of the test or the drug is not present.

Invalid: If a coloured line is not visible in the Control Zone, the test is invalid. Another test should be used to re-evaluate the specimen. Each strip in the Test Cup or Dip is read and functions independently. In the case of an invalid result, the same sample can be used in a second, follow-up test.

An invalid result on one test strip does not invalidate other results derived from the same device.

NOTE: There is no meaning attributed to line colour intensity. Any evidence of a line should be considered a line. In order to prevent incorrect interpretation of the results, the drug test strips should not be interpreted after 10 minutes; adulterant and alcohol results not after 2 minutes.

CUPS

1

Collect specimen into cup

2

Peel the privacy label

3

Read results in 5 minutes.

DIPS

1

Remove cap, dip device into the collected urine for 15 seconds

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