09/12/2017

Tasklist:

- Prepare 0.1M buffer solutions of pH 10, 11, and 12
- Make a 10mM HAuCl₄ and a 100 mM Amino Acid Solutions
- Prepare buffer for Bradford Analysis
- Prepare buffer for Dialysis
- Prepare buffer for SDS-PAGE

My group did not finish the buffer for the dialysis or the SDS-PAGE. These two buffers were prepared the following classes.

▶ 0.1M Buffer Solutions of pH 10, 11, and 12

-Prepare a solution with an acid and its conjugate base by dissolving the acid form of buffer into about 60% of the volume of H2O required to obtain final volume

-Acid: Na2HPO4

-Conjugate Base: Na3PO4

-The objective of the calculations is to find the amount of Na2HPO4 and Na3PO4 in grams to measure out and add into the volumetric flask.

--> Google pKa value: 12.32

--> Equation to use: pH= pKa + log([OH]/[H])

--> [H]= 0.1M since this is a 0.1M buffer solution

--> Calculations for pH 10

Step 1: Find the [OH] pH= pKa + log([OH]/[H]) 10=12.32 + log([OH]/0.1) 0.0047863= [OH]/0.1 [OH]= 4.7863 x 10^-4 M

Step 2: Use [OH] to find the grams of Na3PO4 4.7863 x10^-4M = moles/volume = mole/ 0.25 L moles = 1.1965x10^-4 moles= grams/ MW of Na3PO4 1.1965x10^-4 moles= g/ 163.94 Na3PO4 grams= 0.019616

Step 3: Use [H] to find the grams of Na2HPO4 0.1M= moles/ volume= mole/ 0.25L moles= 0.025 moles= grams/ MW of Na2HPO4 0.025 moles= g/ 141.96 MW Na2HPO4 g= 3.549

--> Calculations for pH 11 and pH 12 were calculated using the same method as above.

> 10mM HAuCl4 and a 100 mM Amino Acid Solutions

--> 10 mM HAuCl4 Solution (Molar Mass of HAuCl4= 339.785 g/mol)

M= mole/ L 10x10^-3M = moles/ 0.25 L Moles= 0.0025

moles= g/MW 0.0025 moles = g/ 339.785 MW HAuCl4 grams= 0.8495

--> Prepare Tubes with HAuCl4, Amino Acid, and Buffer Note: 3 samples were made for each amino acid because the three different buffers (pH 10, pH 11, and pH 12) were used.

	Amino Acid	HAuCl ₄	Buffer
Arg	1 mL	1 mL	8 mL
Phen	1 mL	1 mL	8 mL
Iso	1 mL	1 mL	8 mL

Below are the pictures of these samples.

Figure 1: AuNP and Amino Acids with pH Buffers




Figure 1 shows the different solutions that contained pH buffers. (a) shows the solutions that contained arginine. The sample with a pH of 12 has a purple hue. The other two samples were completely clear. (b) shows the solutions that contained phenylalanine. All three samples were completely clear. (c) shows the solutions that contained isoleucine. All three samples were completely clear.

Remake HAuCl₄ and Arginine Samples with Varying Concentrations of Arginine

	Arginine	HAuCl ₄	Buffer
Test Tube 1	0.01 mL	1 mL	8.9 mL
Test Tube 2	0.1 mL	1 mL	8.9 mL
Test Tube 3	1 mL	1 mL	8 mL

These test tube sets were completed three different times using buffers pH 10, 11, and 12. Below are the pictures of these samples. Unfortunately, we forgot to take the UV/Vis

Figure 2: AuNP and Arg at Different Arg Concentrations and Different pH




Figure 2 shows AuNP and arginine at different arginine concentrations and at different pHs. (a) shows 0.01 mM of arginine at pHs of 10, 11, and 12. The solution at a pH of 11 had a dark purple hue. The solution at a pH of 12 had a light pink hue. (b) shows 0.1 mM of arginine at pHs of 10, 11, and 12. The solution at a pH of 12 had a slight pink hue, and the solution at a pH of 10 had a dark purple hue. (b) shows 2.5 mM of arginine at pHs 10, 11, and 12. All these solutions were completely clear.

Prepare buffer for Bradford Analysis

-1 L 50 mM Tris -50 mM NaCl -pH= 7.5

-Calculation:

(0.05 M)(121.14 MW Tris) = **6.057 g of Tris**

(0.05M(58.44 MW NaCl)= **2.922 g of NaCl**