

Freezing CHO Cells  
Christina McGee 4-5-2012

- 0) Grow 10 cm dish to 50-90% confluency in T-75 flask
- 1) Digest cells with 3 ml Trypsin/EDTA, wait for cells to detach.
- 2) Add 7 ml media, transfer to 15 ml sterile tube.
- 3) Centrifuge 250 G/5 min.
- 4) Remove supernatant. Gently resuspend cell pellet in 2 ml freezing medium\*.
- 5) Transfer 0.5 ml to 4 freezing vials labeled with name of cell line, date and your initials.
- 6) Place vial in freezing box in -80°C freezer.
- 7) Transfer to liquid nitrogen 1-4 days later.
- 8) Record location of cells in Liquid N2 excel spreadsheet

\*Freezing medium: Cell medium (F-12 or DMEM/F-12) without selection agents or antibiotics plus 20% FBS and 10% DMSO (dimethylsulfoxide, sterile).