

# CAMRI MRI MANUAL

Manual Last Updated May 2020

General Information	i
<b>Safety Related Items in the MR Environment</b>	
Metal Screening	1
Ear Plugs and Headphones	2
Medical Crash Cart	2
Subject Squeeze Bulb	3
<b>Safety Policies</b>	
Door Security	5
Accurate Entry of Subject Height, Weight, Age and Sex	5
Temperature Control	5
Obese or large subjects	6
Children	
MRI Principles of MR Imaging and Scanning	7
<b>Safety Related Procedures</b>	
Routine Reboot/ Restart, Full Shutdown and Turning Scanner on procedure	8
Emergency Magnet Quench	11
Emergency Electrical Shutdown	12
<b>Handling Emergencies or Problems</b>	
Medical Emergencies	13
Fire Emergencies & Table Stop	14
Audible Alarms	15
Subject Tingling or Muscle Twitches	16
Subject Perspiration or Increased Heart Rate	17
<b>Special Hazards</b>	
Laser Localizer Light	17
MRI Phantom Leaks	17
Echoplanar (fMRI) Imaging	17
Common unsafe devices chart & Device lookup procedure	18

## GENERAL INFORMATION

### *Risks associated with the MRI lab:*

Used properly, the magnetic resonance imaging equipment contained within the MRI lab is quite safe, however, it poses serious risks to the unwary. Users of the lab should be completely familiar with this manual and with the procedures for protecting others from hazards. To minimize risks to subjects and other members of the research team, **only personnel who have successfully completed the full CAMRI safety certification process are allowed access to the MR scan rooms, control rooms or equipment rooms.** Observers who have not been safety trained are **not permitted** to enter the MRI suites without special prior arrangements.

The main hazards in the lab are:

- The “projectile effect” when heavy, sharp, or dangerous objects are hurled into the instrument. Even seemingly innocuous objects, such as hand tools, can be *lethal*.
- Pacemaker damage. Certain cardiac pacemakers can be damaged by exposure to magnetic fields, causing direct hazards to subjects. Under no circumstances should persons with pacemakers enter the MRI suites at CAMRI.
- As in many laboratories, the MRI lab contains wiring and circuitry that operate at dangerous voltages. Under no circumstances should users touch any exposed wiring, or any exposed terminals in the equipment cabinets.
- Grossly improper scanner operation could result in excessive heating of the subject due to RF energy being deposited. This is easily avoided by operating the equipment according to the guidelines contained in the user manuals and set by the individual instructors.
- Suffocation: in extreme cases, the imaging magnet may release large volumes of helium gas that can rapidly force all air out of the scan room. Normally, the helium gas would be vented through the roof. However, there is a small but significant risk that the venting system could fail.

### *Safe Areas*

There are no areas in the MR suites that can be considered completely safe. The control rooms, scanner rooms, and equipment rooms all have risks associated with magnetic fields and/or electrical equipment. CAMRI safety certification is required for personnel to enter any of these areas.

### *Reduction of Risks*

The chief risk exposure in the lab is to entering personnel who are unfamiliar with the equipment and its hazards. Personnel working in the facility must be constantly vigilant of who is entering the console or scan room areas. Especially in emergency situations, you must ensure that no one without proper training enters any of the scanner rooms, and even then, that they have adequately checked themselves for possible hazards such as projectiles.

Many objects in the scanner control rooms and equipment rooms are **NOT MR compatible** and may become projectiles in the MR scanner rooms. You must never move any object from these rooms into the MR scanner rooms unless you are absolutely certain that the object is MR safe.

Similarly, some objects in the MR scanner rooms may only be safe when kept at a distance from the MR scanner. **Only personnel explicitly authorized to do so should move objects in the scanner room that are labeled “Not MR safe”.** Only objects that are not ferromagnetic should be labeled with a “MR safe” label and this safe label should not be in red or orange. Unlabeled objects should be assumed **NOT** safe to move unless they are clearly non-metallic.

#### *Tours and Training Exercises*

As interesting as the equipment is, please resist the temptation to show visitors the scanner “up close,” as this introduces the unnecessary risk of unwittingly exposing people to potential hazards. Tours or training exercises that would involve having non-safety trained personnel present in the scanners, control rooms or equipment rooms must be authorized in advance by a CAMRI staff member and must be performed in compliance with any special requirements included as part of that authorization.

#### *Reporting of safety incidents or near-incidents*

All incidents or near-incidents must be reported to Krista Runge or Lacey Berry as soon as possible and no more than 24 hours after the incident. Contact information is available at the end of this manual.

Information in this manual that you believe to be incorrect or out-of-date should be reported to Lacey Delay.

In any emergency, try to step through the following guidelines methodically and carefully. Review the safety features, policies and procedures in this manual regularly to assure that you do not need to take unfamiliar actions in a panic situation.

A copy of this manual is available on our wiki: <https://openwetware.org/wiki/CAMRI:Training>

## METAL SCREENING

Anyone preparing to enter an MR scanner room must complete a metal screening form, and this form must be reviewed before access to the scanner room is granted. Individuals who are safety certified at CAMRI are not required to personally complete a formal written metal screening form about themselves but are responsible for verifying that they are personally safe to enter the scanner room.

If there are any doubts regarding the metal screening responses, **do not allow the individual to enter the scanner room**. The fact that the individual has been scanned in an MR scanner previously (even at CAMRI) is **never** a sufficient basis upon which to conclude that the subject can enter the scanner room safely, since risks vary according to magnetic field strength. Dental fillings and orthodontic braces do not constitute significant risks (though the latter may produce unacceptable artifacts) and do not preclude scanning. Subjects with tattoos or permanent eyeliner should be advised of the small risk of local redness or irritation and asked to report any discomfort immediately. Scanning should be stopped immediately if such discomfort develops.

Some implanted metal devices have been established as safe for MR scanning. A recent copy of Sherlock's book cataloguing implanted medical devices is available in the MR suite and up-to-date information is always available on the website <http://www.mrisafety.com>. If, in reviewing these resources, you believe that it is possible to safely scan your subject, you should contact Krista Runge or Lacey DeLay to request authorization to scan the subject. Even if you are certain that the implanted metal does not constitute a risk, **do not allow the individual into the scanner room unless you have obtained explicit authorization to do so**. It is better to ask these questions before your research subject arrives to prevent cancellation at the last minute and inability to scan.

Before entering the scanner room, subjects and staff must remove all objects from their person that might constitute a risk in the MR environment. It is the investigator's responsibility to assure that this has been done. Subjects should be asked to turn pockets inside out to demonstrate that no potentially hazardous objects have been overlooked. Gold and platinum jewelry is not ferromagnetic. Nonetheless, subjects should remove ALL jewelry before going in the scanner since these metals can still conduct electricity and therefore pose a risk for burns in the presence of time-varying magnetic fields. Jewelry that forms large loops is particularly hazardous.

## **EAR PLUGS AND HEADPHONES**

Anyone in the scanner room while the scanner is in operation must be provided with and must use hearing protection in the form of earplugs and/or headphones to avoid hearing injury from the acoustic noise generated by the scanner. We provide both ear plugs and headphones in each of our rooms here at CAMRI.

### Sources of acoustic noise:

When the scanner is not acquiring images, the only source of acoustic noise is the chirping of the compressor pump used to circulate cold water. Different studies have report this noise to be in the 65-80 dB range, below the OSHA standard necessitating the use of a hearing protection device (HPD). Make sure that you always hear this chirping sound before you start scanning.

The time-varying gradient magnetic fields are the primary source of acoustic noise during MR imaging. The rapid rise and fall of currents within the gradient coils in the presence of the static magnetic field create strong Lorentz forces that cause the gradient coils to move against their mountings. The vibration of the coils and the vibration and flexing of their mountings cause the loud tapping and knocking noises during imaging


## **MEDICAL CRASH CART**

We do not have a MD here at CAMRI and we do not have a medical crash cart. For these reasons individuals that pose a risk for coding or any other safety concern should not be scanned here at CAMRI. Because we have no crash cart we do not do studies that involve contrast. In case of an emergency you should dial 911 and remove the individual from the scanner before they arrive. If dialing 911 from scanner please dial a 9 first. (9-911)

## SUBJECT SQUEEZE BULB

### Siemens Squeeze Ball

The scanners are equipped with a squeeze bulb that allows the subject to set off an audible alarm to attract the operator's attention. The squeeze bulb should be made available to subjects unless some alternative method of constant monitoring (e.g., another person in the scanner room) is in effect. Use of the squeeze bulb or some comparable form of continuous subject monitoring is mandatory if you are operating the scanner. The squeeze bulb plugs into the red connector at the foot of the bed. You can verify that the squeeze bulb is connected by noting that the squeeze bulb LED lights up when you press the talk button on the intercom. If the subject squeezes the squeeze bulb, a continuous audible alarm is emitted via the intercom and the intercom squeeze bulb LED lights up.




The silence alarm button for the squeeze ball is labeled "(11)". The Table stop button is labeled "(5)".

(INTERCOM IN THE TRIO SCANNER)



The silence alarm button for the emergency squeeze ball is labeled (1). The emergency table stop button is labeled (2).

(INTERCOM IN THE PRISMA SCANNER)



Squeeze ball connector can be found at the foot of each table. Attach hose to color coded plug if it becomes disconnected.

### **Responding to a squeeze bulb alarm**

1) If a scan is ongoing, press the “Stop” button on the CONSOLE using the mouse. (You will need to do this because the scan is not stopped with the participant squeezes the emergency ball, it only alerts you that it has been squeezed).

*This is especially important if you are scanning in “Level 1” mode rather than standard mode since the subject may have triggered the squeeze ball in response to magnetostimulation or excessive heating.*

2) Press the appropriate intercom talk button to stop the audible alarm. (If you have not stopped the scanner first than you will hear the scanner noises and you will not be able to hear the subject talk)

3) While holding down on the intercom talk button, speak to the subject to determine why the squeeze bulb was pressed. *Make sure that the volume is turned up so that you can hear the subject’s response.*

4) If necessary, enter the room to further investigate and/or correct the problem.

## **POLICY REGARDING DOOR SECURITY**

Scan doors should **NEVER** be left open unattended. Make it a habit to close the door after you. Only level 2 users will have access to the scan room doors. Do not let someone use your badge to access the scanner. ANYONE that is not level 2 is not allowed to be in the scan rooms unattended for any reason, they must be accompanied by a trained individual.

## **POLICY REGARDING ACCURATE ENTRY OF SUBJECT HEIGHT, WEIGHT, AGE AND SEX**

The scanners require that the subject's height, weight, age and sex be entered before scanning.

Accurate information regarding weight must be provided to ensure that FDA limits for energy deposition are not exceeded. Weights should be correct to within five pounds. Incorrect weight should **never** be entered in an effort to get the scanner to conduct a study that it otherwise would not perform because FDA limits would be exceeded. Since we are a research facility we must not use patient identifiers for HIPAA violations and age, sex and height can be anonymized. You should NEVER use the subject's name when registering the subject for a scan at CAMRI. All experiments should have a subject number used in place of the subject's name.

## **POLICY REGARDING TEMPERATURE CONTROL**

In regulating energy deposition in the subject's body in accordance with FDA guidelines, the scanners assume that the ambient temperature in the room does not exceed 72° and that the relative humidity does not exceed 60%. Consequently, the thermostat should never be set for a room temperature higher than 72°. Please note that only cotton, linen or paper should be used for bed covering or blankets since radiofrequency energy may cause heating of synthetic sheets or blankets.



## **POLICY REGARDING OBESE OR LARGE SUBJECTS**

Subjects weighing more than 300 pounds should not be scanned at CAMRI. The Trio 3.0 Tesla scanner bed is designed to support weights up to 440 pounds. Even subjects weighing substantially less than 300 pounds should never be allowed to sit at the distal end of either of the scanner beds, since they are not designed to support the full weight of a large subject applied at full mechanical advantage. The Bore size is 60cm so many patients around 300 lbs may not fit in our bore.

To avoid burns or peripheral nerve stimulation, a minimum distance of 5 mm should be maintained between the subject's body and the wall of the scanner tunnel. MR pads or cotton sheets available in the MR scan rooms can be used to assure this distance is maintained or that skin to bore contact is not being made.

## **POLICY REGARDING CHILDREN**

Children may only enter the scan rooms as participants in an IRB approved research study of children. Children not involved in the research study (e.g, the child or sibling or a research subject) may not enter the scan room and must remain in the waiting room.

All safety precautions applicable to adult subjects are applicable and if anything, more important in children. Careful metal screening, accurate entry of age, sex and weight, and use of "Standard Mode" scanning options whenever possible are important steps in minimizing risks to this population .

### **3 Principles of MR Imaging:**

**THE SCANNER IS ALWAYS ON.** A serious hazard for MRI safety is allowing an object to become attracted to the static magnetic field. This can result in an individual being struck, injured or trapped against the magnet by the object. Equipment can also become damaged by slamming into the magnet or being struck by another object that is accelerating rapidly due to the strong attraction of the magnetic field. There are three main component of the MRI scanner: the static magnetic field, the radiofrequency field, and the gradient or time varying magnetic field.

#### **3.1 Static Magnetic Field**

The static magnetic field is the main magnetic field that is always present once the scanner is ramped up to the designated field strength. This field is described in units of Tesla (T) or Gauss ( $1\text{T} = 10,000 \text{ gauss}$ ). Our system uses a 3T static magnetic field, approximately 60,000 times stronger than the earth's magnetic field that induces a compass to point North. The distance for the magnet that is safe for the general public and to use all objects and devices is denoted as the 5 gauss line. Do not bring anything ferromagnetic or having the ability to become magnetized into the scan room.

#### **3.2 Radiofrequency Field**

The radiofrequency (RF) coil is the heating source within the scanner. This system uses coils as transmitters to excite the MRI signal and as a receiver to detect the MRI signal. It is important to properly use the equipment and accessories of the MRI scanner.

#### **3.3 Gradient/Time varying Magnetic Field**

The gradient or time varying magnetic field selects the slices and imaging planes. This particular field is superimposed over the static magnetic field, and is the sources of all the acoustic noise. The coils within this system are pulsed on and off to produce linear gradients of the magnetic field for imaging. This allows producing an array of images with different spatial and temporal resolutions, and with different contrast between tissues in the image.

### **Performing a Routine Restart or Reboot**

A Restart/ Reboot of the MRI scanner initiates a routine electrical reboot should a situation or problem arise. Your participant can stay in the scanner for this as long as they give you consent. This takes a few minutes to complete and can be done if any of the following occurs:

1. Computer screen is freezing or loading slow.
2. The scanner table is not responding to controls.
3. An error message has occurred that requires the system to be rebooted.

### **Restart or Reboot Procedure: (5-7 Minutes)**

1. Along the top of the system is a bar of various commands. Select “System”.
2. From “System” <End Session>:
3. Then select, <Restart System>
4. Once the scanner beeps 3 times you are able to continue your scan.

Restart System is only used to reboot the system.

### **Performing a Routine Full Shutdown**

A full shutdown of the MRI scanner initiates a routine electrical shutdown should a problem arise. The subject is not allowed to stay in scanner for this and the bed should be brought out to the home position (all the way up and out). The full shutdown process takes approximately 30 minutes and should be done if any of the following occurs:

1. Reboot/Restart did not fix error messages.
2. There are visible alarms that revealed that the magnet has quenched or that helium levels are low.
3. You do not hear the chirping sound or Siemens advises you to turn it off.
4. An error message occurred that requires a full shutdown of the scanner.

### **Full Routine Shutdown Procedure: (30 minutes)**


1. Along the top of the system is a bar of various commands. Select “System”.
2. From “System” <End Session>:
3. Then select <Shutdown System>

This is used when the Restart/ Reboot does not fix error message. If error still persists after full shutdown, contact Siemens.

## PERFORMING A FULL ROUTINE SHUTDOWN

### Routine Shutdown Procedure


1. Click the System tab at the top of the screen:



2. Click on "End Session":



3. Click "Shutdown System".



4. In the confirmation dialog box that appears, click "Yes"
5. It will take approximately 5-10 minutes before you see "It is now safe to turn off your computer" (picture a.) on the MR computer screen. You must then press the blue "system off" button on the Siemens scanner control box (picture b.) located on the wall next to the MR Scanner window.



When it is appropriate to restore power, press the blue “power on” button located on the Siemens control box next to the control room window. (If turning back on after full shutdown please wait 20 seconds before turning system back on) The system will need approximately 20 minutes to reinitialize. To avoid subsequent problems, make sure that the bed is completely out of the gantry at its home position and that the top head coil is attached and plugged in before restoring power. Wait until you hear the 3 loud beeps before moving the table or trying to scan.

#### **To turn on the scanner after shutdown:**

1. Turn the key on control panel to unlock position and press blue scanner on button.

Scanner will take about 20 minutes to boot up before you can use it.

(You will see a red line at bottom of screen while booting)

2. Once you hear 3 beeps then the scanner is ready to use. (The red line will disappear as well) PLEASE NOTE\*\* If you start moving table before 3 beeps the scanner might give you an error message and you will need to perform the 30 minute shutdown to resolve this.

## PERFORMING AN EMERGENCY MAGNET QUENCH

### (LIFE THREATING EMERGENCY involving Magnetic field)

Users of the CAMRI facility should only quench the magnet in the event that the magnetic field itself poses an **immediate risk to an individual's life**. Two such circumstances are:

- 1) A metal object is lodged in the scanner in a way that poses an immediate serious threat to a person (e.g., the person is pinned to the magnet by a metal object that is causing internal injuries).
- 2) Fire personnel determine that there is ***no other alternative*** to entering the room with axes or other heavy gear when fighting a fire.


If the absence of a major emergency, facility users should never quench the magnet by themselves, even if they are convinced that a magnet quench will ultimately be necessary (e.g., if a large object has been drawn into the magnet, but poses no immediate risk to a person). If no one is harmed by a metal object flying into the scanner a controlled quench can be performed by a Siemens engineer that will reserve the helium and allow them to bring down the field in a controlled setting.

### Quench Procedure

1. Locate and press the **QUENCH BUTTON** in the control room or scanner room. The magnetic field will fall to a safe level within 30 seconds.

In the control rooms, quench buttons are located on the top portion of the Siemens wall mounted control boxes located just to the left of the window. Both control room quench buttons are covered by a plexiglass cover labeled “Quench Button For Emergency Use Only”. The quench button itself has the word “STOP” printed three times around its perimeter.

In the scanner room, quench buttons are located behind the door and is represented by the yellow image below.



Trio control room quench button is the button labeled “(1)”



Scan room quench button located behind door

2. When the magnet is quenched, the helium in the scanner boils off. Although the helium should vent out of the room to the rooftop, **you should make sure the door to the scanner**

**room is wide open before quenching the magnet.** If possible, you should remove yourself and the subject from the scanner room before quenching the magnet to minimize the chance of asphyxiation in the event that the helium is improperly vented.

Page 11

3. Dial 9-911 and request medical assistance as detailed elsewhere in this manual.
4. The helium vent ducts become dangerously cold during a quench. Do not touch them.
5. Immediately notify a CAMRI staff member that you have quenched the magnet.
6. At best, it will be at least a few weeks before the scanner can be returned to service. If uninjured, your research subject should be sent home

### **PERFORMING AN EMERGENCY ELECTRICAL SHUTDOWN**

The following events should prompt an emergency electrical shutdown:

1. You see smoke or fire coming from the scanner, equipment room or console.
2. Flooding has carried or is threatening to carry water into electrical equipment

**Electrical shutdowns do not turn off the magnetic field—the magnet is always on unless the magnet has been quenched.**

#### **Emergency Electrical Shutdown Procedure**

1. Locate and press one of the large red electrical shutdown buttons in the scanner room or control room. **Make sure that it is the electrical shutdown button, not the quench button.**



The emergency electrical shutdown button is red and it is labeled Emergency electrical shutdown button.

2. Electrical shutdown immediately stops all power to the scanner, the scanner equipment and the console computers. It does not turn off the lights. Also, power to the stimulation equipment will not be interrupted, so be aware that electrical or fire hazards may still be

present.

3. In the case of fire or medical emergency, call 9-911

#### Page 12

4. Remove your subject from the scanner room. The electric brakes on the scanner bed will have been released automatically, so simply pull the bed out of the gantry manually using the handle at the foot of the bed.
5. Notify CAMRI staff that you have performed an Emergency Electrical shutdown.
6. Circumstances that justify an emergency electrical shutdown do not typically justify quenching the magnet. **Do not quench the magnet unless there is a specific reason to do so** (possible reasons for quenching the magnet are discussed elsewhere in this manual).
7. If uninjured, send your subject home. It will take at least a couple of days to restore the scanner to operational status.

### HANDLING MEDICAL EMERGENCIES

The following procedures are designed on the assumption that a physician or nurse is not immediately available in the MR laboratory at the time of the emergency. If a physician or nurse is present, the medical recommendations may be adjusted as deemed medically appropriate for the subject's condition. However, all non-medical aspects of these guidelines, **particularly those related to removing the person from the magnet or scanner room, must be followed** to avoid unnecessary injury to the subject or personnel.

1. If (**and only if**) the medical emergency involves the subject being pinned to the magnet by a metal object held in place by the magnetic field, quench the magnet following the procedure described elsewhere in this manual.
2. Call 911 (9-911) Describe the event. Advise the person taking the report that the building is a secure building and that you will provide access via the back door, which is the entry closest to you.
3. If the emergency involves a subject in the magnet:
  - A. Stop the scan and go in the room to help subject.
  - B. *Pull the bed completely out of the scanner bore. **The scanner bed cannot be detached from the scanner.***
  - C. Remove the subject from the scanner room so that emergency personnel can safely assess the subject.

***Under no circumstances should a code team or emergency personnel untrained in MR safety enter the scan room. Always remove the subject from the instrument first!***



5. Provide medical assistance in accordance with your training and experience while awaiting arrival of the paramedics. Consider the following options:

A. Initiate CPR if the person is pulseless or not breathing

Page 13

### **HANDLING FIRE EMERGENCIES**

1. In case of fire, call 911 (9-911 if calling from CAMRI phone).
2. If smoke or fire is coming from the scanner, equipment room or console, perform an emergency electrical shutdown as described elsewhere in this manual. You will need to manually pull the bed out of the scanner because power to the scanner will be lost. To do this hit the table stop button (either on scanner or on intercom), this will release the bed from the tracks and allow you to manually pull bed out.
3. If you are scanning and smoke or fire is NOT coming from the scanner, equipment room or console, stop the scan and safely bring your person out of the scanner. If time permits, initiate a routine electrical shutdown by selecting “End” from the “System” menu at the far right at the console.



The Table STOP button is labeled “(1)”

**WE DO NOT HAVE MRI SAFE FIRE EXTINGUISHERS HERE AT CAMRI**

4. Remove the subject from the scanner and escort the subject out of the building.
5. Do not return to the building until advised by fire personnel that it is safe to do so.
6. Contact CAMRI personnel to advise them that there was a fire in the building.

### **TABLE STOP BUTTON**

A table stop is performed when you need the bed to stop because something is being pulled or stuck in the bed tracks. It will release the bed from the guide tracks and allow you manually pull

the bed out. The table stop can also be used if you need to manually pull someone out of the scanner if you performed an electrical shutdown, bed is not responding or if power is off.

## Releasing the table stop

Scanner 5 (TRIO) –

From the Scanner: Go into the room and press the Table up and table down buttons in succession. This will reset the table and allow the bed to be used again.

From the console: Click on table properties button and click release table stop button. This will reset the table and allow the bed to be used again.

Scanner 4 (PRISMA)-

From Scanner: Go into the room and turn the table stop button (on bed) clockwise and then press the table up and down buttons in succession. This will reset the table and allow the bed to be used again.

From the console: Hit the small flashing button on back of intercom and go into room and press the table up and table down buttons in succession. This will reset the table and allow the bed to be used again.

## HANDLING AUDIBLE ALARMS


You should never scan while an audible scanner-related alarm is sounding. If you cannot identify and correct the underlying problem, your study should be discontinued. If an audible alarm is sounding, investigate the following possibilities:

1) The alarm might be the building fire alarm. This extremely loud alarm is audible throughout the building, is associated with flashing lights in the hallways, and would be difficult to mistake for a scanner related alarm. Even if you suspect that the fire alarm has been triggered accidentally, you **MUST** do the following:

- A. If you are scanning, stop the scan.
- B. Go into the scanner room and pull the bed out.
- C. Assist the subject off the bed
- D. Accompany the subject out of the building via the nearest accessible exit.
- E. Do not reenter the building until told that it is safe to do so by fire personnel.

2) The alarm might have been triggered by someone squeezing the squeeze bulb. Look to see if the squeeze bulb LED on the intercom is lit. If it is, see the separate section regarding the squeeze bulb. You will be able to continue your study if this is the source of the alarm.

3) The helium level might be low or the magnet might have quenched spontaneously due to an earthquake or as a result of someone pressing the quench button. Check the Siemens control box located in the control room immediately to the right of the window. If the magnet stop LED is lit, the magnet has quenched. If the helium level LED is lit, the helium level is low. On the Trio, other LEDs indicate power problems, compressor problems, battery problems and communications errors. You can press the alarm silencer to stop the audible alarm, but **do not scan**. Notify CAMRI staff of the problem and send your subject home.



Trio

The magnet stop LED is labeled “(3)”

The helium level LED is labeled “(2)”

The alarm silenced (acknowledge) button is labeled “(9)”

4) The alarm might be a building related alarm. Check the annunciator panel in the hallway near the building exit for an LED indicating the source of the problem. If the LED's indicate a

problem that is outside the MR suite, you will generally be able to continue your study. Facilities should automatically be notified of alarms appearing on this board.

Page 15


### **Audible Noise Check:**

The scanner makes a **chirping sound**, this sound is perfectly normal and should be heard on all properly functioning scanners while not in use. This chirping sound is produced by the compressor and if you do not hear this sound please let CAMRI staff know. If the sound is not present it could be the chiller not functioning properly and the scanner may heat up. The scanner will be out of commission for repairs if the compressor is down.

## **HANDLING SUBJECT TINGLING OR MUSCLE TWITCHES**

Tingling or muscle twitches are potential physiologic effects of time varying magnetic fields. Such effects are particularly likely to occur with echo-planar imaging in fMRI studies. To minimize the likelihood of such magnetostimulation, operate the scanner in “Standard Mode”. In this mode, only 1% of subjects should experience such effects. However, the scanner may refuse to scan certain subjects with certain pulse sequences in “Standard Mode”. If you operate in “Level 1” operating mode, up to 50% of subjects may experience magnetostimulation with certain pulse sequences.

Complaints of tingling or muscle twitches should prompt rescreening for any metal objects that might have been previously overlooked and verification that subject positioning does not form potential loops. For echo planar imaging, selecting a phase encoding direction that is anterior-posterior (when this is an option) should reduce the likelihood of magnetostimulation. Note that the sensory input associated with magnetostimulation will pose an unwanted confound in fMRI studies.



Subject positioning loops that predispose to magnetostimulation or burns

### **HANDLING PERSPIRATION AND INCREASED PULSE AND SUBJECTS WITH CONDITIONS ASSOCIATED WITH IMPAIRED THERMAL REGULATION**

Perspiration and an increased pulse rate may result from energy deposition in the body during scanning. Energy deposition in the body is carefully regulated by the scanner in accordance with FDA guidelines. If your subject develops these symptoms, you should verify that the subject's age, height and weight were entered correctly when registering the patient, since these parameters may influence the calculated energy deposition. You should also verify that the room temperature does not exceed 72° and the humidity does not exceed 60% since the calculated energy deposition limits assume that they do not. The Trio will measure the temperature and may refuse to scan certain sequences if the temperature exceeds 71.6°. For subjects who have medical conditions such as fever, diabetes, pregnancy, or cardiovascular disease that can impair thermal regulation, you should operate the scanner in "Standard Mode" if possible, since energy deposition is not a concern in this mode. Children or elderly subjects are also at increased risk of overheating. If you do scan subjects with conditions associated with impaired thermal regulation in "Level 1" mode, you should be attentive to signs or symptoms of overheating and stop the scan if overheating is suspected. "Level 1" mode should be avoided if possible in subjects who are unable to communicate reliably (e.g., children, sedated subjects, stroke patients). Adjusting the fan in the scanner may be helpful in reducing the likelihood of overheating in subjects.

### **LASER LIGHT LOCALIZER HAZARDS**

On the 3.0 Tesla Trio scanners, a laser is available for landmarking the patient's position in the scanner. Subjects should be instructed to keep their eyes closed while the laser light is turned on to avoid eye injury. If the laser light appears as a spot rather than as crosshairs, it should be turned off immediately, and you should notify one of the designated CAMRI staff that it is in need of repair.

### **MRI PHANTOM LEAK HAZARDS**

The MR phantoms used to calibrate the scanners are sealed and should not show any evidence of leakage. The contents of some of the phantoms is potentially hazardous. If a phantom develops a leak, protective clothing (gloves, labcoat, goggles and, if the contents have become aerosolized, a face mask) should be worn while cleaning the leak. The contents should be disposed of as hazardous materials (i.e., not simply poured down the drain).

### **ECHOPLANAR (fMRI) IMAGING HAZARDS**

Echoplanar imaging, used in fMRI, utilizes rapidly changing gradients and is associated with higher voltages than many other MR imaging modalities. The risk of magnetostimulation is increased with echoplanar imaging. The risk of magnetostimulation can be reduced by choosing a phase-encoding direction that is oriented anterior-posterior when this is an option.

### The List of Common Unsafe Devices

• Cardiac Pacemaker	• Vascular access port and/or catheter*
• Implanted cardioverter defibrillator (ICD)	• Radiation seeds/implants
• Electronic Implant or Device	• Swan-Ganz (thermodilution catheter)
• Magnetically-activated implant/device	• Mediation patch (Niction, Nitroglycerin)
• Neurostimulator system (Deep Brain)	• Any metallic fragment/foreign body
• Neurostimulator system (Urinary Control)	• Wire mesh implant
• Spinal cord stimulator	• Tissue expander (e.g. breast)
• Internal electrodes/wire	• Surgical staples, clips, or metallic sutures*
• Bone growth	• Joint replacement (hip, knee, etc.)*
• Cochlear, otologic, or ear implant	• IUD, diaphragm, or pessary
• Insulin or infusion Pump	• Dentures or partial plates*
• Implanted drug infusion pump	• Tattoo or permanent makeup
• Prosthesis (eye and/or penile)	• Body piercing jewelry
• Heart valve prosthesis	• Hearing aid*
• Eyelid spring/wire	• Other implant
• Artificial or prosthetic limb	• Breathing problems or motion disorder
• Metallic stent, filter, coil*	• Claustrophobia
• Shunt	• Aneurysm Clips

\*\*\* This is not a complete list of unsafe objects, just the most common\*\*\*

### **Contra-indicators of MR**

There are equipment, devices, medical implants, etc that are at risk to injury and/or damage if they are entered into the scanner room. This pertains to be within the body or in personal possession. The items that are listed are considered to be potential risk to an individual's safety. Some items on the list are conditional (\*), and further information can be founded at [www.mrisafety.com](http://www.mrisafety.com)

