


# Dexcom CLARITY<sup>®</sup>

Dexcom CLARITY<sup>®</sup> software captures your CGM data so you can view glucose patterns, trends and statistics anytime, anywhere. Upload glucose data from your Dexcom CGM receiver or auto-stream with the G5 Mobile system. Share data reports with your clinic in ways that work best for you.



## Access Dexcom CLARITY on the Internet

- 1** Enter **clarity.dexcom.eu** in your internet browser.
- 2** Click **Dexcom CLARITY<sup>®</sup> for Home Users**.
- 3** Enter your username and password for any of these accounts:  
Dexcom Account, Dexcom G5<sup>®</sup> Mobile App or Dexcom CLARITY
- 4** If you do not have any Dexcom account:  
Click **Create Account** and complete the registration form.  
Click **Sign in**.

While logged into CLARITY, click **Help** at the top of the page to read the User Guide and learn more about CLARITY features.

For help and support with Dexcom CLARITY® software, call Technical Support at 1-844-832-1809.

Review CLARITY computer and internet requirements at [www.dexcom.com/clarity-requirements](http://www.dexcom.com/clarity-requirements).



# Dexcom®



Dexcom, Inc.  
6340 Sequence Drive  
San Diego, CA 92121 USA  
+1 858.200.0200  
[www.dexcom.com](http://www.dexcom.com)

Dexcom Canada, Inc.  
501 - 4445 Lougheed Hwy.  
Burnaby, BC V5C 0E4 Canada  
Sales: 1-844-832-1810  
Technical Support: 1-844-832-1809

#### Intended Use/Safety Statement

The web-based Dexcom CLARITY® software is intended for use by both home users and healthcare professionals to assist people with diabetes in the review, analysis, and evaluation of historical CGM data to support effective diabetes management. It is intended for use as an accessory to Dexcom CGM devices with data interface capabilities. Caution: The software does not provide any medical advice and should not be used for that purpose. Home users must consult a healthcare professional before making any medical interpretation and therapy adjustments from the information in the software. Caution: Healthcare professionals should use information in the software in conjunction with other clinical information available to them.

Dexcom, Dexcom CLARITY® and the Dexcom G5® Mobile CGM System are either registered trademarks or trademarks of Dexcom, Inc. in the United States and/or other countries. All other product or company names that may be mentioned in this publication are tradenames, trademarks or registered trademarks of their respective owners.

LBL015367 Rev001