

- Announcements
- Quiz (last one!)
- Pre-lab Lecture
 - ❖ DMMB assay
 - ❖ More about ELISA
 - ❖ qPCR analysis
 - ❖ Today in Lab (M3D6)

Announcements

- Lec 7: Atissa on proposal; Lec 8: special topics
- Next time: finish report + clean up; WAC survey
- Mod 3 research proposal
 - feedback on your wiki page today/tonight
- Mod 3 report
 - no separate Methods section needed: *do* state any unique conditions along the way in Results section/captions
 - required analysis: viability (incl. stats), qPCR (incl. RNA amount and cell counts), ELISA (incl. standard curves)
 - optional/if relevant analysis: PG assay; general bead, cell, and media appearance on D2, D3, and/or D4

DMMB assay

- Measure GAGs with cationic dye
- Absorbance shift due to complex $A_{595} \downarrow$
 - fades quickly! (pipet upstairs)
 - at low pH, selects for sulfates over carboxyls
 - thus a correction for alginate
 - standard curve made with chondroitin-6 sulfate
- Typically normalize to cell amount (cf DNA content in qPCR)... maybe next year


ELISA protocol

- Direct ELISA uses labeled primary antibody
- Indirect ELISA – why use a secondary antibody?
 - Signal amplification
 - flexibility/efficiency (t, δ) → use w/many 1° Abs
 - but may be more cross-reactive
- Development process – what/why/how
 - 2° Ab – AP conjugated
 - provide substrate: pNPP → colorimetric rxn.
in soln; A₄₂₀
 - * development time is key
 - detect low [protein] → long enough
 - ↓ not oversaturate everything → not too long

ELISA Outcomes

Outcome	Possible Explanations
High reading in "blank" samples	<ul style="list-style-type: none">• cross-contamination of wells or global in PB S• stopped & lock step• incomplete washes
No signal at all (including standards)	<ul style="list-style-type: none">• wrong or old/bad Abs• flipped plates at 1st Abs step• too high [Tween]
Saturated signal for some samples <u>Egypt'1</u>	<p>too concentrated → run dilution series usually</p>

ELISA (and PG) analysis



qPCR analysis

efficiency of each primer set \times

$$\text{ratio} = \frac{(E_{\text{target}})}{(E_{\text{ref}})}^{\Delta CP_{\text{target}}(\text{control-sample}) - \Delta CP_{\text{ref}}(\text{control-sample})}$$

$E_{\text{target}} = 2.0$ $\Delta CP_{\text{target}} = 2$ ratio = $\frac{2^2}{2^0} = 4$ -fold change
 $E_{\text{ref}} = 2.0$ $\Delta CP_{\text{ref}} = 0$

* we can know expression of CNII A vs. CNII B
 BUT NOT CNII A vs. CNI A
 in absolute sense

Equation 1 from M.W. Pfaffl, Nucleic Acids Res 29:2002 (2001)

Today in Lab (M3D6)

- Finish ELISA – includes 90 min incubation
- Meanwhile... *figure out how much DEV SOLN your super-group needs!*
 - DMMB assay staggered (15-20 min of work)
 - qPCR analysis
 - finish viability analysis if you haven't already
 - cross-group research discussion
 - (optional: start clean-up)