


- Announcements
- Lab Quiz
- Pre-lab Lecture
 - ❖ Plasmids: Basics + Gel EP
 - ❖ DNA Ligation
 - ❖ Bacterial Transformation
 - ❖ Today in Lab: M2D4

Announcements

- Module 1 drafts returned tonight
 - 2 wks to revise, due 11 am Thu/Fri
- Module 2 versus Module 1
 - growing sophistication and self-reliance
 - in lab work *protocols will be less explicit*
 - in report *focus on core results (not each step)*
 - this doesn't mean don't ask questions!

lac operon

These three genes encode metabolic enzymes



Encodes a repressor protein that binds to P_{LAC} (operator) turning it OFF.

I_{ATG}
↑
In turn, if lactose binds to the repressor, it is made inactive, turning ON expression of P_{LAC}.

Plasmid Overview

antibiotic resistance




origin dependent
affect copy #

↳ circular, ds, extrachromosomal

why? vector for inserting
a gene into an organism

ampicillin resistance
→ select bacteria that
have plasmid w/ Amp

DNA EP: Shape-dependence



Plasmid versus linear samples

e.g., 4 kbp plasmid

1-cut: linear, runs w/ ladder


2-cut: sums to 4 kbp

uncut: supercoiled - faster

circular, (relaxed, nicked): s lower

+ high MW concatemers

Restriction Enzymes for Cloning



CCGGCC = G
ccCTAG

pED-IPTG-INS
oligo → ins
7 bkb


What if BamHI is 5' and XmaI is 3'
on insert? **insert + reversed,
non-coding product**

What if BamHI on 5' and 3' ends?
non-directional cloning, bkb re-close

Can you get multiple inserts?

-X X B X X B B-

DNA Ligation



Reaction creates **phosphodiester bond**

Reaction requires **ATP**

What factors affect yield?

ratio of substrates

[DNA], rxn. time


enzyme (conc., T, age, conc.)

How do we assess if it worked?

sequencing

diagnostic digest if possible


Bacterial transformation



Other methods: electroporation, ballistics

DNA Amplification in Bacteria

Plate bacteria + DNA mixture



Transformation Controls + Outcomes

NExt TIME

| Sample | Expectation... What if? | Role |
|--------|----------------------------|------|
| | | |
| | | |
| | | |
| | | |

Miller assay calculation and result

$$1 \text{ Miller Unit} = 1000 * \frac{(Abs420 - (1.75 * Abs550))}{min - t * v * Abs600}$$

x2 if plate assay
if stock OD600, use stock v01.
!! diluted

What was the key finding from your first assay?

What do you expect from your new construct?

Today in Lab (M2D4)

- Decide on a workflow you like
 - Goal A: Ligation/transformation of pED-IPTG-YFD
 - Goal B: IPTG → lacZ transfer function (plate assay)
- FNT
 - I will send an email when it's finalized
 - Have a great spring break!