- Announcements
- Pre-lab Lecture
  - Writing a methods section
  - Where we are/going (cloning)
  - DNA cleanup
  - Today in Lab: M1D2

### **Announcements**

- Brief discussion of orientation day quiz ★C V = CQ V2 \*
- General reminders

- Assignments + Schedule = our syllabus 
$$\% = \frac{3}{2} / \frac{1000}{1000}$$

- All lecture notes linked from Schedule page, "notes"
- Assignments (scroll down) has notebook guidelines
- First M1 quiz next time! (see wiki)
- · Optional ungraded notebook submission today "PLEASE GRADE" @ for MIDZ entry by Sat SIM

# Mod 1 major assessments

- Methods section (5%)
  - Done individually, no revision
  - Important building block/skill for Mod 2 report
  - Plasmid construction only
  - Start in today's FNT!
- Abstract and data summary (15%)
  - Done in pairs, with revision
  - More building blocks: abstract, figures/captions, etc.
  - But not yet a full narrative
  - Starts in next FNT and today's WAC session


# Methods section tips

- Organizing sub-sections
   in a log: row order
  - Often start with an overview/introductory sentence (what, not why) → then give step-by-step details
- Methods should be <u>concise</u> and complete at a protocol
  - Space-wise, avoid tables/lists when a sentence will do
  - Sentence-wise, avoid extra words
  - Content-wise, cover what's needed and only that needed to understand and replicate your work
- Concentrations are more useful than volumes; or you can state amounts, plus total volume.


### Methods section exercise

- Consider the following passage: "Template DNA (5 ng) and primers were mixed with 20 uL of 2.5X Master Mix in a PCR tube. Water was added to 50 uL. A tube without template was prepared and labeled control."
  - What information is missing?
  - What information can be cut?


#### Roadmap for Plasmid Construction




### Where we are





# Digesting PCR product






## Restriction enzymes for cloning





### Restriction enzyme analysis example



Enzyme	# bands	Size 6
EcoRI	2	~700 ~4800
EcoRI + Xbal	3	~ 4700
EcoRI + BamHI	4	~ 300 doublet

Salt +EtOH GRH waste no bleach!



Silica resin column

[qiagen.com]

DNA clean-up

1. bind DNA TSulf, UpH chartropic salt disrupts H-bonds, delydrates DNA, bridges to silica

2. keep DNA, woshe Ise EtoH keeps precipitated

3. elute DNA Jeat, high pH
repulsion -Si 10- 5-1-0

## Today in Lab: M1D2

- use 10ml (not 4)
- Careful with enzyme stocks!
  - Keep cold; don't contaminate
- There are 4 samples today (2 single digests)

  PCR-dbl NNY-dbl XbI Ecorl-HF (+wcu+)
- Order of addition for digest
- WAC session at 3:45 pm
- Note: FNT uses ApE useful general DNA tool todder quoted "Do you have eyeloshes in there?

  Science proposal: Con I take teen?" NOT IRB-approved