


System Engineering
 20.109(F09)
 M2D3 lecture
 10.22.09


Roadmap for System Engineering Module


Day	Lab	Lecture
✓ 1	Testing v1.0	System Eng: Bact. Photography
✓ 2	Measuring	Two Component Signaling
3	Tools	Synthetic Biology + Genetics
4	Journal Club	Group Meeting
5	Re-tune	Examining the C-dog
6	Readout v2.0	System Eng: revisited
7	Journal Club	Group Meeting

New tools for reliable engineering of complex biological systems




Synthetic Biology


Standardization 


Abstraction 


Synthesis A+A+C+T+T...


Genetic Engineering

rDNA 

Sequencing 

PCR 





Registry of Standard Biological Parts

Part:BBa_115010
 Designed by Jeff Tabor Group: UTAustin (2004-09-20)

cph8 (Cph1/EnvZ fusion)
 Chimeric Cph1 light receptor/EnvZ protein.


Usage and Biology
 Requires PCB biosynthetic genes (BBa_115008 and BBa_115009) for chromophore formation. The light responsive domain (Cph1) has maximal response to light near 660nm. Exposure to red light inhibits the activity of the EnvZ histidine kinase domain. When active (in the dark) EnvZ phosphorylates endogenous OmpR, a transcription factor which activates transcription from the OmpC promoter and represses transcription from the OmpF promoter. Must be used in E.coli deficient in native EnvZ.


New tools for reliable engineering of complex biological systems

Does the Registry have Cph8 K⁺ or P⁺ signaling mutants?

Looking for parts: genetic screening

Step 1: Mutagenize gene of interest
 Step 2: Put DNA in cells (if not there already)
 Step 3: Look for mutant phenotype





Library Variations

Cph1/EnvZ	A553	G554	V555	S556	H557
EnvZ	A239T	G240E	V241G	S242D	H243A
wt seq	GCG	GGG	GTA	AGT	CAC
oligo seq	RNS	RNS	RNS	RNS	SNW
poss aa	Val	Val	Val	Val	Val
	Ala	Ala	Ala	Ala	Ala
blue = K-P+	Asp	Asp	Asp	Asp	Asp
red = K+ P-	Glu	Glu	Glu	Glu	Glu
	Gly	Gly	Gly	Gly	Gly
	Ile	Ile	Ile	Ile	Leu
	Met	Met	Met	Met	Pro
	Thr	Thr	Thr	Thr	His
	Asn	Asn	Asn	Asn	Gln
	Lys	Lys	Lys	Lys	Arg
	Ser	Ser	Ser	Ser	
	Arg	Arg	Arg	Arg	

NOTE: no stop codons should be in mix

