



RIVERBAND[®]

SURGICAL ADHESIVE



PORTLAND, OR USA
www.rpmed.com

RiverBand™ Surgical Adhesive

Premium Veterinary Glue



New 3mL EasySqueeze™ Bottle.

RiverBand™ bottle is redesigned in a larger size to hold more glue and provide more dispensing comfort.



RiverBand™
3mL Size



Vibrant Blue Color Identification.

RiverBand™ contains purified n-Butyl Cyanoacrylate with a D&C certified dye that allows for differentiation between the product and bodily fluid. Surgical Adhesive should be applied topically; it is not intended for internal closures.



RiverBand™
Vibrant Blue



Competitors
Dull Purple



Sterile & Accurate Applicators.

The disposable applicators come in a sterile pack, which compared with the competition, can make a difference in an effort to reduce bacterial migration into the surgery site. The tip fits snugly into the tip on the bottle and allows greater control over the size and placement to the drop of product.



RiverBand™
Sterile Pouch



Competitors
Non-Sterile Pouch



Fast Dry Time.

RiverBand™ Surgical Adhesive dries clear for clean appearance. It sets up quickly (5-15 seconds) and dries flexible for maintaining closure on wound site.



RiverBand™
5+ seconds



Competitors
8+ seconds

RiverBand™ | Offered now in a 3mL bottle!



RiverBand™

Designed for wound management
in situations such as:

- ✓ Lacerations and abrasions
- ✓ Declaws
- ✓ Cuts
- ✓ Tail docking




Order Codes

| Item Number | Description | Size |
|-------------|-------------------|------|
| TA-03 | Surgical Adhesive | 3 mL |

Package includes:

- 3.0 mL dropper bottle
- 12 sterile disposable applicator tips
- Instructional insert





CORPORATE HEADQUARTERS

825 NE 25th Avenue
Portland, OR 97232 USA

OFFICE +01 (503) 517-8001

TOLL FREE +01 (866) 445-4923

FAX +01 (503) 517-8002

WEB www.rpmed.com | **EMAIL** info@rpmed.com

ISO 13485:2003 | ISO 9001:2008 | FDA Registered

**MANUFACTURING IN PORTLAND,
OREGON SINCE 1990.**

Riverpoint Medical designs, develops and markets next generation medical devices for the Veterinary and Healthcare markets.

Specialized product lines include Surgical Sutures, MedLED® Cordless Surgical Lighting, Oncology Products, and Orthopedic Sutures and Devices.