



MedLED[®]

LIGHTING



PORTLAND, OR USA
www.rpmed.com

Welcome to the MedLED Chrome®

The next frontier in the elimination of cords in the O.R.

Finally, the gap is bridged between Xenon power and LED efficiency. The **Chrome** is the first & brightest LED headlight that provides portable power all day long.

The **Chrome** utilizes the flagship *truly cordless* **MedLED®** headstrap design for complete freedom. By choosing among a family of removable light modules, the **Chrome** can be tailored to perform any surgery, even those extremely demanding deep cavity cases where brightness is imperative.



The 4 MedLED® Pillars

1. Be free of the cords with the patented battery-swap technology giving extensive battery life.



2. A patented headband system to hold everything together in perfect balance. So much, that you forget it's there.



med
LED[®] **CHROME**[®]
TRULY CORDLESS™



3. At 275,000 LUX, the **Chrome**[®] MC7 PRO is the brightest surgical LED headlight in the world.



4. Tailor your light. Choose between a soft-top or hard-top headband and load on an MC3, MC6, or MC7 PRO.

The MedLED Chrome® Family



MC3 Chrome® (new Onyx®)

E E-CLASS *Clinical*

- Dental & Hygiene
- Perio Clinic
- Oral & Maxillofacial Clinic
- Plastic & Cosmetic Clinic
- OB/GYN Clinic
- ENT Clinic
- Veterinary Clinic



MC6 Chrome®

C CROSSOVER

- Plastic, Cosmetic &
- Oral & Maxillofacial
- Neuro Surgery
- Vascular Surgery
- OB/GYN Surgery
- ENT Surgery
- General Surgery
- Veterinary Surgery



(new Sapphire[®])

Clinical & O.R.

Recon. Surgery
Surgery

MC7 PRO Chrome[®]

P PREMIUM *O.R.*

- Cardio & Cardio-Thoracic Surgery
- Plastic, Cosmetic & Recon. Surgery
- Spine Surgery
- Oral & Maxillofacial Surgery
- Ortho Trauma Surgery
- Orthopedic Surgery
- Neuro Surgery
- General Surgery
- Urology

MC3 Chrome® (new Onyx®)

E E-CLASS *Clinical*



MC6 Chrome® (new Sapphire®)

C CROSSOVER *Clinical & O.R.*



MC7 PRO Chrome®

P PREMIUM *O.R.*




Chrome® Specifications vs Others*




SPOT SIZE

Metric


MC3



MC6




**MC7
PRO**




BATTERY LIFE

Hours



BRIGHTNESS

LUX



*data on file

MedLED®

A Product Line by Riverpoint Medical
Manufactured in Portland, Oregon USA



TRULY CORDLESS™
Premium Surgical Headlights

MedLED® Story *The headlamp that started a revolution.*

MedLED® was born as the first in surgical LED headlights and remains as the only *truly cordless™* option for any surgical application.

It all started in 2004 when we ‘saw the light’ and launched the **MedLED® Classic** - a bright, soft strap LED headlight followed by the **MedLED Focus®**, **Onyx®** and **Sapphire®**. **And now, the MedLED Chrome®.**

Headlights

MC3	100,000 LUX
MC6	200,000 LUX
MC7 PRO	275,000 LUX

Headstraps

Hard-Top (HT) Top Coronal Adjustment Knob
Soft-Top (ST) Soft Neoprene Velcro Top

Kit Contents

Standard Kit (SK) Light, 2x Batt., USB/Plug & Case
Hospital Kit (HK) Light, 6x Batt., 2x USB/Plug, 4-Bay Charger, +2D Magnifiers & Case.

MC3 Chrome® (new Onyx®)



MC3-HT-SK	Hard-Top	Standard Kit
MC3-ST-SK	Soft-Top	Standard Kit
MC3-HT-HK	Hard-Top	Hospital Kit
MC3-ST-HK	Soft-Top	Hospital Kit

MC6 Chrome® (new Sapphire®)



MC6-HT-SK	Hard-Top	Standard Kit
MC6-ST-SK	Soft-Top	Standard Kit
MC6-HT-HK	Hard-Top	Hospital Kit
MC6-ST-HK	Soft-Top	Hospital Kit

MC7 PRO Chrome®




MC7-HT-SK	Hard-Top	Standard Kit
MC7-ST-SK	Soft-Top	Standard Kit
MC7-HT-HK	Hard-Top	Hospital Kit
MC7-ST-HK	Soft-Top	Hospital Kit

MC3/MC6/MC7 PRO Chrome®




MC367-HT-SK	Hard-Top	Standard Kit
MC367-ST-SK	Soft-Top	Standard Kit
MC367-HT-HK	Hard-Top	Hospital Kit
MC367-ST-HK	Soft-Top	Hospital Kit




Onyx®
Sapphire®


Chrome®




Wall
Plug



4-Bay
Desktop Charger



Onyx® / Sapphire® / MC3
Small Clip Magnifier



Chrome® MC6 / MC7 PRO
Large Clip Magnifier



Magnifier
Example

Accessory

ORDER CODE	DESCRIPTION	WORKS WITH
------------	-------------	------------

Batteries

M-90002	Battery Set (2 pcs.)	Onyx / Sapphire
M-90004	Battery Set (2 pcs.)	Chrome

Chargers


M-90102	Wall Plug & USB Cable (US)	All Units
M-90106	Wall Plug & USB Cable (US/EU/UK/AU)	All Units
M-90151	USB Cable (Only)	All Units
M-90202	4-Bay Desktop Charger & USB Cable (US)	All Units
M-90203	4-Bay Desktop Charger & USB Cable (US/EU/UK/AU)	All Units

Magnification


Diopter (Magnification, Working Distance)

M-90301	+2.0D (40% Increase, 250mm Distance)	Onyx / Sapphire / MC3
M-90302	+2.5D (50% Increase, 235mm Distance)	Onyx / Sapphire / MC3
M-90303	+3.0D (60% Increase, 220mm Distance)	Onyx / Sapphire / MC3
M-90304	+3.5D (70% Increase, 210mm Distance)	Onyx / Sapphire / MC3
M-90305	+4.0D (80% Increase, 200mm Distance)	Onyx / Sapphire / MC3
M-90306	+4.5D (90% Increase, 195mm Distance)	Onyx / Sapphire / MC3
M-90321	+2.0D (40% Increase, 250mm Distance)	MC6 / MC7 PRO
M-90322	+2.5D (50% Increase, 235mm Distance)	MC6 / MC7 PRO
M-90323	+3.0D (60% Increase, 220mm Distance)	MC6 / MC7 PRO
M-90324	+3.5D (70% Increase, 210mm Distance)	MC6 / MC7 PRO
M-90325	+4.0D (80% Increase, 200mm Distance)	MC6 / MC7 PRO
M-90326	+4.5D (90% Increase, 195mm Distance)	MC6 / MC7 PRO


Padding



Snap Padding




Pillow Velcro (Hard-Top)




Pillow Velcro (Soft-Top)

M-90901	Neoprene Snap Set (3 pcs.)	Onyx / Sapphire
M-91302	Neoprene Snap Set (3 pcs.)	Chrome Hard-Top / Soft-Top
M-91305	Pillow Velcro Set (2 pcs.)	Chrome Soft-Top
M-91306	Pillow Velcro Set (3 pcs.)	Chrome Hard-Top
M-91312	Neoprene Snap 10 Sets (3 pcs. ea.)	Chrome Hard-Top / Soft-Top
M-91315	Pillow Velcro 10 Sets (2 pcs. ea.)	Chrome Soft-Top
M-91316	Pillow Velcro 10 Sets (3 pcs. ea.)	Chrome Hard-Top


Other




Carrying Case



GoPro Mount




Amber Filter




Headband Insert

M-90912	Carrying Case	Chrome
M-90921	Camera Mount for GoPro® Hero 3 & 4 Black	Chrome
M-90922	Camera Mount for GoPro® Hero 5 Session	Chrome
M-90932	Camera Mount for LoupeCAM® HD	Chrome
M-91101	Amber UV Filter for Composite Curing Prevention	Onyx / MC3
M-91103	Amber UV Filter for Composite Curing Prevention	MC6
M-91331	Headband Insert Velcro for Additional Frame Support	Chrome

Headstraps & Bezels



Headstrap Hard-Top



MC3 MC6 MC7 PRO

M-92301	Headstrap Hard-Top (HT)	Chrome
M-92302	Headstrap Soft-Top (ST)	Chrome
M-92331	Bezel MC3	Chrome
M-92361	Bezel MC6	Chrome
M-92371	Bezel MC7 PRO	Chrome

Headstrap Soft-Top Example



Hospital Kit Example





CORPORATE HEADQUARTERS

825 NE 25th Avenue
Portland, OR 97232 USA

OFFICE +01 (503) 517-8001
TOLL FREE +01 (866) 445-4923
FAX +01 (503) 517-8002
WEB www.rpmed.com
EMAIL info@rpmed.com

ISO 13485:2003
ISO 9001:2008
FDA Registered

PUB304-03 12-2017

**MANUFACTURING IN PORTLAND,
OREGON SINCE 1990.**

Riverpoint Medical designs, develops and markets next generation medical devices for the Veterinary and Healthcare markets.

Specialized product lines include Surgical Sutures, MedLED® Cordless Surgical Lighting, Oncology Products, and Orthopedic Sutures and Devices.



**RIVERPOINT
MEDICAL**