



NEEDLE TYPES

REVERSE CUTTING ▾

3/8 Circle




1/2 Circle




TAPERSHARP™ ▾

3/8 Circle



1/2 Circle



PREMIUM COSMETIC CONVENTIONAL CUTTING ▲

3/8 Circle




CONVENTIONAL CUTTING ▲

1/2 Circle



PREMIUM REVERSE CUTTING ▾

3/8 Circle



STRAIGHT CUTTING ▾











KS
60mm















NEEDLE TYPES

TAPER POINT

3/8 Circle

- BV50-2
2.6mm 
- BV70-3
3.8mm 
- BV100-3
3.8mm 
- BV5-100-4
4.7mm 
- BV100-5
5.1mm 
- BV140-5
5.1mm 
- BV140-6
6.4mm 
- BV-130-5
6.5mm 
- BV-1
9.3mm 
- BV175-8
9.3mm 

1/2 Circle

- BV130-5
6.5mm 
- TF
13mm 
- RB-1
17mm 
- SH-1
22mm 
- CT-3
22mm 
- SH
26mm 
- M0-6
26mm 
- CT-2
26mm 
- CT-1
36mm 
- CT
40mm 
- CTX
48mm 
- TP-1
65mm 



BLUNT POINT

1/2 Circle

- CTB
40mm 

SPATULA


1/4 Circle

- S-14
8mm 
- S-29
8mm 

3/8 Circle

- TG160-6
5.5mm 
- TG140-6
6.5mm 
- TG140-8
6.5mm 
- CS140-6
6.5mm 

1/2 Circle

- CS160-6
5.5mm 
- CS160-8
5.5mm 