



ONCOLOGY PRODUCTS



PORTLAND, OR USA
www.rpmed.com

BRACHYTHERAPY NEEDLES

ROUND HUB NEEDLES

SQUARE HUB NEEDLES



ROUND HUB




SQUARE HUB

ADVANTAGES

- Electropolished for sharpness
- Siliconized for easy insertion
- Mirror polished I.D. for a smooth, easy seed deployment
- Echo-enhanced tipped for ultrasound visualization
- Centimeter markings for depth placement
- Rounded stylet for better seed delivery
- Foil hub marker to indicate bevel tip orientation

| PRODUCT CODE | HUB | DESCRIPTION | FOIL MARKER COLOR | STYLET HUB COLOR |
|---------------|--------|--------------------|-------------------|-------------------|
| RP-1120-1820 | ROUND | 18g (1.2mm X 20cm) | Blue foil marker | Blue stylet hub |
| RP-1100-1820 | SQUARE | 18g (1.2mm X 20cm) | Blue foil marker | Blue stylet hub |
| RP-1100Y-1820 | SQUARE | 18g (1.2mm X 20cm) | Blue foil marker | Yellow stylet hub |
| RP-1100-1812 | SQUARE | 18g (1.2mm X 12cm) | Blue foil marker | Blue stylet hub |

APPLICATOR NEEDLES



APPLICATOR

ADVANTAGES

- Electropolished for sharpness and siliconized for easy insertion
- Mirror polished I.D. for a smooth, easy seed deployment
- Echo-enhanced tipped for ultrasound visualization
- Centimeter markings for depth placement
- Trocar stylet for smooth direct entry
- Designed to fit the Mick® 200 T.P.V. Applicator

| PRODUCT CODE | DESCRIPTION |
|--------------|--------------------|
| RP-1200-1820 | 18g (1.2mm X 20cm) |
| RP-1200-1720 | 17g (1.4mm X 20cm) |

BRACHYTHERAPY ACCESSORIES

BRACHYTHERAPY GRID



ADVANTAGES

- Easy Needle Insertion
- Precise 18g Needle Channel Guides
- Accurate Placement
- Single, One-Time Sterile Use
- Ensure that Ultrasound Tip Matches Grid Graphics

| PRODUCT CODE | DESCRIPTION |
|--------------|-------------|
| RPG18 | 18 Gauge |

BIOPSY NEEDLES

MAGNUM STYLE

BIOPTY / PRO-MAG*** STYLE



ADVANTAGES

- Etched Cutting notch for better tissue retention
- Echo enhanced tip for ultrasound placement
- Spacer clip for easy loading
- Centimeter markings for depth placement
- Electropolished for sharpness & smooth entry
- 19mm or 17mm cutting notch for tissue capture

| PRODUCT CODE | STYLE | HUB | DESCRIPTION | FOIL MARKER COLOR | STYLET HUB COLOR |
|--------------|------------------|-------|--------------------|-------------------|------------------|
| RP-1000-1820 | MAGNUM | ROUND | 18g (1.2mm X 20cm) | Blue foil marker | Blue stylet hub |
| RP-1000-1825 | MAGNUM | ROUND | 18g (1.2mm X 25cm) | Blue foil marker | Blue stylet hub |
| RP-1120-1820 | BIOPTY / PRO-MAG | ROUND | 18g (1.2mm X 20cm) | Blue foil marker | Blue stylet hub |

GOLD FIDUCIAL MARKERS & KITS

**BIOPTY® is a registered trademark of BARD®
***PRO-MAG™ is a trademark of ANGIOTECH

SQUARE HUB GOLD MARKER KIT



ADVANTAGES

- Electropolished for sharpness
- Siliconized for easy insertion
- Mirror polished I.D. for a smooth, easy seed deployment
- Echo enhanced tipped for ultrasound visualization
- Centimeter markings for depth placement
- Rounded stylet for smooth gold marker delivery
- Foil marker to indicate bevel tip orientation

| PRODUCT CODE | DESCRIPTION |
|--------------|--|
| GM 1231-1720 | (1) 1.2 X 3mm Gold Marker with 17g x 20cm Needle |
| GM 1233-1720 | (3) 1.2 x 3mm Gold Marker with 17g x 20cm Needle |
| GM 1251-1720 | (1) 1.2 X 5mm Gold Marker with 17g x 20cm Needle |
| GM 1253-1720 | (3) 1.2 x 5mm Gold Marker with 17g x 20cm Needle |
| GM 1031-1820 | (1) 1.0 X 3mm Gold Marker with 18g x 20cm Needle |

| PRODUCT CODE | DESCRIPTION |
|---------------|--|
| GM 1033-1820 | (3) 1.0 x 3mm Gold Marker with 18g x 20cm Needle |
| GM 1051-1820 | (1) 1.0 X 5mm Gold Marker with 18g x 20cm Needle |
| GM 1053-1820 | (3) 1.0 x 5mm Gold Marker with 18g x 20cm Needle |
| GM 0831-1820* | (1) 0.8 x 3mm Gold Marker with 18g x 20cm Needle |
| GM 0833-1820* | (3) 0.8 x 3mm Gold Marker with 18g x 20cm Needle |

COUPLED GOLD MARKER KIT



| PRODUCT CODE | DESCRIPTION |
|--------------------|--|
| GM 0831-SS17-1820* | 1) 0.8 x 3mm Gold Marker with 18g x 20cm Needle |
| GM 0832-SS17-1820* | (2) 0.8 x 3mm Gold Marker with 18g x 20cm Needle |
| GM 0831-SS20-1820* | (1) 0.8 x 3mm Gold Marker with 18g x 20cm Needle |
| GM 0832-SS20-1820* | (2) 0.8 x 3mm Gold Marker with 18g x 20cm Needle |

[Additional combinations may be available on request; Gold markers also available seperately]

*0.8 markers are not knurled



CORPORATE HEADQUARTERS

825 NE 25th Avenue
Portland, OR 97232 USA

OFFICE +01 (503) 517-8001

TOLL FREE +01 (866) 445-4923

FAX +01 (503) 517-8002

WEB www.rpmed.com | **EMAIL** info@rpmed.com

ISO 13485:2003 | ISO 9001:2008 | FDA Registered

**MANUFACTURING IN PORTLAND,
OREGON SINCE 1990.**

Riverpoint Medical designs, develops and markets next generation medical devices for the Veterinary and Healthcare markets.

Specialized product lines include Surgical Sutures, MedLED® Cordless Surgical Lighting, Oncology Products, and Orthopedic Sutures and Devices.