

Sanitizer Test Log

Registro de Sanitizante

Southwest Training Systems
520.744.1092 / 1.888.838.1550

School: _____ Manager: _____

- Immerse Test Papers in STILL solution for 10 seconds
(Sumerja Papeles de Prueba en solución por 10 segundos)
- Buckets & Sinks must be @ 200 ppm *minimum* (200 to 400 is acceptable)
(Cubetas y Fregaderos deben medir 200ppm. 200 ha 400 es acceptable)
- Bottles must be @ 200 ppm *minimum* (Bottellas deben medir 200ppm. 200 ha 400 es acceptable)
- Strength complies with Maricopa Health Code Chap 8, Regulation 8, Section R (3)(e)(1.4) – 21CFR 178.1010 (b)22 & (c)17 - 40CFR 153-168-180

Example
Ejemplo

Sanitation Procedures

Date	Time	Source Test	ppm Test	Action	Signature
01/01/01	12:00	  	200 +	N/A	Mr. Clean
		  			
		  			
		  			
		  			
		  			
		  			
		  			
		  			
		  			
		  			
		  			
		  			
		  			
		  			
		  			



Bucket



Bottle



Sanitize Sink